

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KLF4 Anti-Human Antibody (Clone AT4E6)

单克隆抗体

产品货号: MA01463

第三版

免疫原:Recombinant human KLF4 (1-170aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:AT4E6

描述:KLF4, which is also known as the gut-enriched Kruppel-like factor (GKLF), is a zinc fingercontaining transcription factor that belongs to the Kruppel-like family of transcription factors. KLF4 plays important roles during the proliferation and differentiation of epithelial cells.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF