



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KRT14 Antibody (Clone AT13F4)

单克隆抗体

产品货号: MA01467

第三版

免疫原:Recombinant human KRT14 (1-472aa) purified from E. coli.

亚型:Mouse IgG2a, k

克隆:AT13F4

描述:Keratin 14 (KRT14) is a member of the type I keratin family of intermediate filament proteins. KRT14 gene provides instructions for making a protein called keratin 14. Keratins are a group of tough, fibrous proteins that form the structural framework of certain cells, particularly cells that make up the skin, hair, and nails. Keratin 14 is specifically produced in cells called keratinocytes in the outer layer of the skin.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF