

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

KRT18 Anti-Human Antibody (AT2G3)

单克隆抗体

产品货号: MA01468

第三版

免疫原:Recombinant human KRT18 (79-430aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:AT2G3

描述:Keratins (cytokeratins) are intermediate filament proteins that are mainly expressed in epithelial cells. Keratin heterodimers composed of an acidic keratin (or type I keratin, keratins 9 to 23) and a basic keratin (or type II keratin, keratins 1 to 8) assemble to form filaments.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF