

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Lag (Langerhans Cells), (4D12) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01474

第三版

免疫原: Human Langerhans Cells

亚型: IgG1

克隆: 4D12

描述: Langerhans cells (LC) were reported to function as the antigen-presenting cells in the skin. They express immune response gene-associated antigens on their cell membranes.

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF