



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LAIR-1 Anti-Human Antibody (Clone lc12)

单克隆抗体

产品货号: MA01475

第三版

免疫原: Recombinant human LAIR-1 (22-125aa) purified from E. coli

亚型: Mouse IgG1 heavy chains and light chain

克隆: lc12

描述: Reacts with extracellular domain of the LAIR-1, a member of the leukocyte-associated Ig like receptor family. The LAIR-1 is a transmembrane glycoprotein (~32kDa) expressed on the majority of human peripheral blood mononuclear leukocytes including T cells, B cells, NK cells, macrophages and dendritic cells

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF