

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LDHA Anti-Human Antibody (AT1A4)

单克隆抗体

产品货号: MA01477

第三版

免疫原: Recombinant human LDHA (1-332aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: AT1A4

描述: Lactate through lactate dehydrogenase A (LDHA), which is encoded by a target gene of c-Myc and hypoxia-inducible factor (HIF-1). LDHA expression is repressed by SugR in the absence of sugar.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF