

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

LEC26 (LA102) Anti-Mouse Rat IgG MoAb

单克隆抗体

产品货号: MA01479

第三版

免疫原: Mouse lymphangioma

亚型: IgG2b,

克隆: LA102

描述: Recent years, LYVE-1 (1999, Banerji et al.), podoplanin (1999, Breiteneder-Geleff et al.) and Prox-1 (1999, Wifflé & Oliver) were found one after another as specific antigens of lymphatic endothelium and the approach to lymphatic study has progressed dramatically. However, it has been reported that these antigens are not expressed in only lymphatic endothelium according to tissues or clinical conditions, but antibodies against these antigens also react with other vessels especially vascular endothelium.

配方: Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF