



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LITAF Anti-Human Antibody (Clone AT5C10)

单克隆抗体

产品货号: MA01482

第三版

**免疫原:** Recombinant human LITAF (1-161aa) purified from E. coli

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** AT5C10

**描述:** Lipopolysaccharide-induced tumor necrosis factor-alpha factor (LITAF) mediates the expression of inflammatory cytokines such as TNF-alpha in Lipopolysaccharide-induced processes. LITAF binds to STAT6B, a member of the STAT6 family forming a complex on the TNF-alpha promoter that modulates TNF activity.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF