



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LTA Anti-Human Antibody (AT15A3)

单克隆抗体

产品货号: MA01486

第三版

免疫原: Recombinant human LTA (35-205aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: AT15A3

描述: Lymphotoxin-alpha (LT-alpha; LTA) is member of the TNF ligand superfamily, bind the same TNF receptor and mediate similar pleiotropic effect. LTA, located within the MHC III region of chromosome 6, show close linkage to the HLA class I (HLA-B) and class II (HLA-DR) genes.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF