

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Lysyl-tRNA synthetase (KRS)

单克隆抗体

产品货号: MA01487

第三版

**免疫原:** Full length native human recombinant KRS.

**亚型:** IgG

**克隆:**

**描述:** Lysyl-tRNA synthetase(KRS) is a component of multi-ARS complex and catalyzes the ligation of lysine with lysyl-tRNA KRS is secreted from cancer cell by TNF-alpha and activates.

**配方:** Anti-lysyl-tRNA synthetase antibody containing 0.05% sodium azide for preservative.

**浓度:** Please see datasheet.

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF