
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

LYVE-1 Anti-Human Antibody (Clone 4G1)

单克隆抗体

产品货号: MA01488

第三版

免疫原: Recombinant human LYVE-1 (25-235 aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: 4G1

描述: Lymphatic vascular endothelial hyaluronan receptor-1 (LYVE-1), a selective marker of the lymphatic endothelium, is a surface endocytic receptor for hyaluronan (HA), which is an extracellular glycosaminoglycan involved in cell adhesion and migration.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF