



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MAPK1 Anti-Human Antibody (AT1A4)

单克隆抗体

产品货号: MA01492

第三版

免疫原: Recombinant human MAPK1 (1-360aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: AT1A4

描述: Mitogen-activated protein kinases (MAPKs) are reported to be critical regulatory factors for the growth and migration of various cell types including vascular smooth muscle cells (VSMCs). It has also been reported that the activation of the MAP kinase family, extracellular signal regulated kinases 1/2 (ERK1/2), and c- Jun N-terminal kinase (JNK) via Src activation is important for AII-induced migration of VSMCs.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF