



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

MAPK3 Anti-Human Antibody (AT1A2)

单克隆抗体

产品货号: MA01494

第三版

**免疫原:**Recombinant human MAPK3 protein (1-137aa) purified from E. coli

**亚型:**Mouse IgG2b heavy chain and kappa light chain

**克隆:**AT1A2

**描述:**Mitogen-activated protein kinase (MAPKs) are a widely conserved family of serine/threonine kinase involved in many cellular programs such as cell proliferation, differentiation, motility, and death. MAPK3 (ERK1) signaling pathway can be activated in response to a diverse range of extracellular stimuli including mitogens, growth factors, and cytokines and is an important target in the diagnosis and treatment of cancer.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在-20 °C 或者-70 °C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF