

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Maspin Anti-Human Antibody (Clone AT2N6)

单克隆抗体

产品货号: MA01496

---

第三版

**免疫原:**Recombinant AGXT (330-392aa) purified from E. coli

**亚型:**Mouse IgG2b heavy chain and kappa light chain

**克隆:**AT2T4

**描述:**The AGXT gene provides instructions for making a liver enzyme called alanine-glyoxylate aminotransferase. Inside liver cells, this enzyme is found in peroxisomes, structures that contain many different enzymes used to produce energy and the basic materials important for cellular activities.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在-20 °C或者-70 °C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF