



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MAT2A Anti-Human Antibody (Clone AT3A2)

单克隆抗体

产品货号: MA01498

第三版

免疫原: Recombinant human MAT2A (1-395aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: AT3A2

描述: Methionine adenosyltransferase II (MAT II) is a key enzyme in cellular metabolism and catalyzes the formation of S-adenosylmethionine (SAMe) from L-methionine and ATP. MAT2A is expressed in extrahepatic tissues. In liver, MAT2A expression associates with growth, dedifferentiation, and cancer.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF