

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

PRDX1 Anti-Human Antibody (Clone AT1D8)

单克隆抗体

产品货号: MA01499

第三版

免疫原:Recombinant peroxiredoxin-1 (1-199aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:AT1D8

描述:PRDX1 is a member of the peroxiredoxin family of antioxidant enzymes, which reduce hydrogen peroxide and alkyl hydroperoxides. This protein is also known to be important to protect red blood cells against reactive oxygen species and in tumor prevention.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF