

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Psma Anti-Human Antibody (Clone K1H7)

单克隆抗体

产品货号: MA01506

第三版

**免疫原:**Recombinant human PSMA (117-351aa) purified from E. coli

**亚型:**Mouse IgG1 heavy chain and kappa light chain

**克隆:**k1H7

**描述:**Prostate-specific membrane antigen (PSMA) is a 100 – 120 kDa type II transmembrane glycoprotein that is a glutamate-preferring carboxy-peptidase. Three functionally distinct proteins are encoded, including folylpoly-gamma-glutamate carboxypeptidase in the intestine, N-acetylated alpha-linked acidic dipeptidase 1 in the brain, and prostate-specific membrane antigen in the prostate

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF