
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MB Anti-Human Antibody (Clone AT6E10)

单克隆抗体

产品货号: MA01507

第三版

免疫原: Recombinant human MB(1-154aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT6E10

描述: MB (Myoglobin) is an iron- and oxygen-binding protein. It is found in abundance in the muscle tissue of vertebrates in general and in almost all mammals. MB (myoglobin) is related to hemoglobin, which consists of four myoglobin-like subunits that form a tetramer and are responsible for carrying oxygen in blood.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF