



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MBD3 (Methyl-CpG-Binding Protein 3) (3A3) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01509

第三版

免疫原:Synthetic peptide of part of C terminal of human MBD3 (EELARDGEAPLDKA)

亚型:IgG1

克隆:3A3

描述:Methylation at CpG dinucleotides in genomic DNA is a fundamental epigenetic mechanism of gene expression control in vertebrates. Proteins with a methyl-CpG-binding domain (MBD) can bind to single methylated CpGs and ,ost of them are involved in transcription control.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF