

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PSMD11 Anti-Human Antibody (Clone AT1F4)

单克隆抗体

产品货号: MA01510

第三版

免疫原:Recombinant human PSMD11 (1-422aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:AT1F4

描述:PSMD11 (Proteasome 26S non-ATPase subunit 11), also known as S9, Rpn6 or p44.5, has multisubunit protease complexes consisting of 20S subunits composed of four seven-numbered rings with two outer rings containing alpha subunits and two central rings composed of beta subunits, and 19S caps of 6 ATPase and 11 non-ATPase subunits.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF