

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-MBP Antibody (Clone 10C12)

单克隆抗体

产品货号: MA01514

第三版

免疫原: Recombinant MBP purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 10C12

描述: Epitope tagging offers an easy and universal strategy for the identification and purification of proteins derived by recombinant DNA technology. The insertion of a Maltose Binding Protein (MBP) tag creates a stable fusion product that does not interfere with the bioactivity of the protein or with the biodistribution of the MBP tagged product.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF