



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MCP-1 Anti-Rat Mouse IgG MoAb

单克隆抗体

产品货号: MA01516

第三版

免疫原:Rat MCP-1 Recombinant Protein

亚型:IgG1

克隆:MB10

描述:Human macrophage/monocyte chemotactic and activating factor (MCAF) also called human macrophage chemoattractant protein-1 (MCP-1) is an 8.6 kDa protein containing 76 amino acid residues. It plays an important role in the inflammatory response of blood monocytes and tissue macrophages.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF