

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

PTS Anti-Human Antibody (Clone AT2B2)

单克隆抗体

产品货号: MA01517

第三版

免疫原: Recombinant human PTS (1-145aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: AT2B2

描述: PTS or PTPS (6-Pyruvoyl tetrahydroprotein synthase) is involved in the biosynthesis of tetrahydrobiopterin (BH₄), the cofactor of several enzymes such as the nitric oxide synthase and the aromatic amino acid hydroxylases.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF