



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RAB27A Antibody (Clone AT2E2)

单克隆抗体

产品货号: MA01521

第三版

**免疫原:** Recombinant human RAB27A (1-221aa) purified from E. coli

**亚型:** Mouse IgG1, k

**克隆:** Anti-human RAB27A mAb, clone AT2E2, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human RAB27A protein

**描述:** Ras-related protein Rab-27A (RAB27A) plays a role in cytotoxic granule exocytosis in lymphocytes. Required for both granule maturation and granule docking and priming at the immunologic synapse.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF