



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MIF Anti-Human Antibody (Clone 4E4)

单克隆抗体

产品货号: MA01524

第三版

免疫原: Recombinant human MIF (1-114 aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 4E4.

描述: The cytokine Macrophage migration inhibitory factor (MIF) has been identified to be secreted by the pituitary gland and the monocyte/macrophage and to play an important role in endotoxic shock. MIF has the unique property of being released from macrophages and T cells in response to physiological concentrations of glucocorticoids.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF