



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MINCLE Anti-Human Antibody (AT16E3)

单克隆抗体

产品货号: MA01525

第三版

免疫原:Recombinant human MINCLE (41-219aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:AT16E3.

描述:MINCLE (macrophage inducible C-type lectin, also called CLEC4E or CLECSF9), which is a diverse family of protein which was originally defined by their ability to recognize a wide range of ligand of carbohydrate structure.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF