



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Rad51 (1B1) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01527

第三版

免疫原:Recombinant Human Rad51

亚型:IgG1

克隆:1B1

描述:The highly conserved *Saccharomyces cerevisiae* Rad51 protein plays a central role in both mitotic and meiotic homologous DNA recombination. Seven members of the Rad51 family have been identified in vertebrate cells, including Rad51, Dmc1, and five Rad51-related proteins referred to as Rad51 paralogs, which share 20 to 30% sequence identity with Rad51.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF