

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RALA Antibody (Clone AT41D5)

单克隆抗体

产品货号: MA01528

第三版

免疫原: Recombinant human RALA (1-203aa) purified from E. coli

亚型: IgG2b,k

克隆: Anti-human RALA mAb, clone AT41D5, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human RALA protein.

描述: RALA (Ras-related protein Ral-A) belongs to the small GTPase superfamily, Ras family of proteins. RALA mediates a distinct downstream signaling pathway from Ras that facilitates cellular transformation. RALA is also thought to be involved in numerous signaling cascades, including regulation of the cytoskeleton, vesicle trafficking, and endocytosis

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF