

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

RAP2B Antibody (Clone AT62D9)

单克隆抗体

产品货号: MA01531

第三版

**免疫原:**Recombinant human RAP2B (1-183aa) purified from E. coli

**亚型:**IgG1,k

**克隆:**Anti-human RAP2B mAb, clone AT62D9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human Rap2b protein.

**描述:**RAP2B is 183 amino acid protein. It belongs to the Ras-related protein family. RAP2B is small GTP-binding protein which cycles between a GDP-bound inactive and a GTP-bound active form. RAP2B is involved in EGFR and CHRM3 signaling pathways through stimulation of PLCE1. It may be associated with cytoskeletal rearrangements and regulate cell spreading through activation of the effector TNIK and regulate membrane vesiculation in red blood cells.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF