

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-mTOR (N5D11) Mouse IgG MoAb

单克隆抗体

产品货号: MA01533

第三版

免疫原:Recombinant Rat mTOR

亚型:IgG1

克隆:N5D11

描述:Rapamycin is a lipophilic macrolide compound and induces a G1 phase growth arrest in *S. cerevisiae* and mammalian cells. The mammalian target of rapamycin (mTOR) has a molecular weight of 290kDa and a Ser/Thr protein kinase activity.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF