

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

MYL2 Anti-Human Antibody (Clone AT3B2)

单克隆抗体

产品货号: MA01538

第三版

免疫原:Recombinant human MYL2 (1-166aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:AT3B2

描述:Myosin, light chain 2 (MYL2) encodes the regulatory light chain associated with cardiac myosin beta heavy chain. It is an important protein involved in the regulation of myosin ATPase activity in smooth muscle and Ca⁺ triggers the phosphorylation of regulatory light chain that in turn triggers contraction.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF