



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RASSF1A Anti-Human Antibody (Clone 3F3)

单克隆抗体

产品货号: MA01539

第三版

免疫原: Recombinant human RASSF1A (1-340aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: 3F3

描述: The tumor suppressor gene RASSF1A (3p21.3), human Ras association domain family 1A, is expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney, pancreas, spleen, thymus, prostate, testis, small intestine, colon, and peripheral blood leukocytes.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF