

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

MYL4 Antibody (Clone AT4E8)

单克隆抗体

产品货号: MA01540

第三版

**免疫原:**Recombinant human MYL4 (1-197aa) purified from E. coli

**亚型:**IgG2a ,k

**克隆:**Anti-human MYL4 mAb, clone AT4E8, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human MYL4 protein.

**描述:**Myosin light chain 4 is a protein that in humans is encoded by the MYL4 gene. Myosin is a hexameric ATPase cellular motor protein. It is composed of two myosin heavy chains, two nonphosphorylatable myosin alkali light chains, and two phosphorylatable myosin regulatory light chains. This gene encodes a myosin alkali light chain that is found in embryonic muscle and adult atria. Two alternatively spliced transcript variants encoding the same protein have been found for this gene.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF