

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RBP4 Anti-Human Antibody (Clone AT2B4)

单克隆抗体

产品货号: MA01541

第三版

**免疫原:** Recombinant human RBP4 (19-201aa) purified from E. coli

**亚型:** Mouse IgG1 heavy chain and kappa light chain

**克隆:** AT2B4

**描述:** Retinol binding protein 4 (RBP4) belongs to the lipocalin family and is the specific carrier for retinol in the blood. This protein was found to be expressed and secreted by adipose tissue, and was strongly associated with insulin resistance.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF