



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

/Repeat Domain (2B11) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01542

第三版

免疫原: Human Tau/Repeat Domain (PGGGKVQIVYKP)

亚型: IgG1

克隆: 2B11

描述: Tau protein is a microtubule-associated protein expressed in neuroaxon and has molecular weight of 50~70 kDa. Some Alzheimer's patients have sites where tau protein is accumulated abnormally.

配方: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度: Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF