



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Nanog Anti-Human Antibody (Clone 5A10)

单克隆抗体

产品货号: MA01546

第三版

免疫原: Recombinant human Nanog (1-154aa) purified from E. coli

亚型: Mouse IgG2a heavy chain and kappa light chain

克隆: 5A10.

描述: Nanog is a multidomain homeobox transcription factor that functions to maintain the undifferentiated state of pluripotent stem cells. Nanog expression counteracts the differentiation-promoting signals induced by the extrinsic factors LIF, Stat3 and BMP

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF