

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

NANS Anti-Human Antibody (Clone AT1G6)

单克隆抗体

产品货号: MA01547

---

第三版

**免疫原:** Recombinant human NANS (1-359aa) purified from E. coli

**亚型:** Mouse IgG 1 heavy chain and light chain

**克隆:** AT1G6

**描述:** NANS(Sialic acids synthase) is an enzyme that functions in the biosynthetic pathways of sialic acids. Expressed ubiquitously, NANS enzymatically catalyzes the H<sub>2</sub>O-dependent formation of Nacetylneuraminic acid (Neu5Ac) and 2-keto-3-deoxy-D-glycero-D-galacto-nononic acid (KDN), both of which are sialic acids.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF