



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Rexo1 Anti-Human Antibody (Clone 4H5)

单克隆抗体

产品货号: MA01548

第三版

免疫原: Recombinant human REXO1 (1060-1221aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: 4H5

描述: REXO1 is an exonuclease domain-containing protein that can bind to Elongin. The Elongin complex stimulates the rate of transcription elongation by RNA polymerase II by suppressing the transient pausing of the polymerase at many sites along the DNA template.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF