
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Napsin A (TMU-Ad02) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01549

第三版

免疫原:Synthetic peptide for a part of N terminus of Human Napsin A

亚型:IgG1

克隆:TMU-Ad02

描述:Recently, a new protein isolated by two-dimensional immuno-electrophoresis from tissue of lung adenocarcinoma was identified as Napsin A, a new aspartic Proteinase. The reactivity of the antibody showed immunohistochemically in type II pneumocytes, alveolar macrophages, renal tubules and exocrine glands and ducts in pancreas.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF