



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RGMa (410-10) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01552

第三版

**免疫原:** Synthetic peptide of a part of Human RGMa (LYERTRDLPGRAAAGL)

**亚型:** IgG2b

**克隆:** 410-10

**描述:** RGMa (Repulsive Guidance Molecule a) is a GIP-anchored protein, about 33 kDa, and it is implicated in axonal guidance and neural tube closure of embryonic retina (ref. 1, 2, 3).

**配方:** Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN<sub>3</sub>

**浓度:** Please see datasheet.

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF