



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

NDE1 Anti-Human Antibody (Clone AT1F7)

单克隆抗体

产品货号: MA01556

第三版

免疫原: Recombinant human NDE1 (1-335aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT1F7

描述: NDE1 is important to microtubule reorganization, which is required for centrosome duplication and the formation of the mitotic spindle and cell migration. NDE1 is a binding partner of DISC1 along with LIS1 at the centrosome and at synapses in neurons.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF