

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

NELL2 Anti-Human Antibody (Clone AT13E7)

单克隆抗体

产品货号: MA01558

第三版

**免疫原:**Recombinant human NELL2 (30-258aa) purified from E. coli

**亚型:**Mouse IgG2b heavy chain and light chain

**克隆:**AT13E7

**描述:**Protein kinase C-binding protein, NELL2 is expressed abundantly in neural tissues and may play important roles in the development of neural tissues. Also, NELL protein contains six epidermal growth factor (EGF)-like repeat domains that participate in calcium binding.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在-20 °C或者-70 °C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF