

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

RNH1 Anti-Human Antibody (Clone AT1H23)

单克隆抗体

产品货号: MA01559

第三版

免疫原:Recombinant Ribonuclease inhibitor 1(RNH1) (7-461aa) purified from E. coli

亚型:Mouse IgG2a heavy chain and kappa light chain

克隆:AT1H23

描述:Ribonuclease inhibitor 1(RNH1) the placental ribonuclease inhibitor, is an acid protein with a molecular weight of 50 kDa. RNH1 is a member of a family of proteinaceous cytoplasmic RNase inhibitors that are expressed in many tissues and bind to both intracellular and extracellular RNases in the cytosol.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF