



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-S-100 Protein (2A10) Mouse Ascites MoAb

单克隆抗体

产品货号: MA01561

第三版

免疫原:Bovine Brain S-100 Protein

亚型:IgG1

克隆:2A10

描述:S100 proteins (16 members) show a very divergent pattern of cell- and tissue-specific expression, of subcellular localizations and relocations, of post-translational modifications, and of affinities for  $\text{Ca}^{2+}$ ,  $\text{Zn}^{2+}$ , and  $\text{Cu}^{2+}$ , consistent with their pleiotropic intra- and extracellular functions.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF