
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Neurofilament, 160K, 210K, (N21-8-4) Mouse Ascites MoAb

单克隆抗体

产品货号: MA01564

第三版

免疫原:Bovine Brain Neurofilament

亚型:IgG1

克隆:N-21-8-4

描述:Neurofilament is intermediate filament which expressed in nerve cells. And the protein consists of three subunits, 200 kDa, 160 kDa and 6.8 kDa. These subunits exists in phosphorylated form in axonal transport. This antibody reacts with the 200 kDa and 160 kDa subunits, and further the antibody is useful for research of neurodegenerative disease.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF