

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Nfatc1 Anti-Human Antibody (Clone At1C3)

单克隆抗体

产品货号: MA01566

第三版

**免疫原:**Recombinant human NFATc1 (428-716aa) purified from E. coli

**亚型:**Mouse IgG2a heavy chain and kappa light chain

**克隆:**AT1C3.

**描述:**Nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 1(NFATc1) is a component of the nuclear factor of activated T cells DNA-binding transcription complex which consists of at least two components: a preexisting cytosolic component that translocates to the nucleus upon T cell receptor (TCR) stimulation, and an inducible nuclear component.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF