

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

SAR1B Antibody (Clone AT1C7)

单克隆抗体

产品货号: MA01567

第三版

免疫原: Recombinant human SAR1B (1-198aa) purified from E. coli

亚型: IgG1,k

克隆: Anti-human SAR1B mAb, clone AT1C7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human SAR1B protein.

描述: SAR1B (GTP-binding protein SAR1b) is 198 amino acid protein. SAR1B is involved in transport from the endoplasmic reticulum to the Golgi apparatus. This protein is activated by the guanine nucleotide exchange factor PREB. It is involved in the selection of the protein cargo and the assembly of the COPII coat complex. Mutations in the SAR1B gene are associated with CMRD which is autosomal recessive disorder of severe fat malabsorption.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF