



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Niban (110C4) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01569

第三版

免疫原:Recombinant human Niban protein

亚型:IgG1

克隆:110C4

描述:Niban has been cloned from cell lines derived from Eker rat by AFLP-Differential Display method as a novel gene which is increased in the early stage of renal carcinogenesis (ref. 6). Then, it has been reported that Niban protein is expressed in various cancer tissues including thyroid cancer. In research of recent years, it has been reported that Niban is an ER stress-inducible protein and it regulates phosphorylation of eIF2 and S6K1/4E-BP1 (ref. 3).

配方:Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF