

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Nicastrin Anti-Human Antibody (Clone AT3B2)

单克隆抗体

产品货号: MA01571

第三版

免疫原: Recombinant human Nicastrin (34-153aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: AT3B2.

描述: Nicastrin is a large protein that is a component of an enzyme called gamma-secretase, which is lodged in the cell's membrane. Nicastrin protein is physically associated with presenilin and plays an important role in stabilization and correct localization of presenilin to the membrane-bound gamma-secretase complex.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF