

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

NKp44 Monoclonal Antibody (Clone AT1G6)

单克隆抗体

产品货号: MA01572

第三版

免疫原: Recombinant human NKp44 (19-130aa) purified from E. coli

亚型: IgG2a, k

克隆: Anti-human NKp44 mAb, clone AT1G6, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human NKp44 protein.

描述: NKp44 (NCR2) is a member of the natural cytotoxicity receptor (NCR) family that is expressed on activated human NK cells. NKp44 displays a single extracellular Ig-like V domain and a transmembrane portion containing the charged residue (Lysine), likely involved in the association with KARAP/DAP12 molecules.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF